



# WILL YOUR LAUNDRY PASS ACCREDITATION?

Australian Laundry Standard AS 4146:2000

This standard requires your linen be disinfected either thermally or chemically.

## THERMAL DISINFECTION

To disinfect linen thermally the wash load must be maintained at a minimum 65°C for 10 minutes.

This requires a main wash cycle time of 14-18 minutes to allow for mixing and to allow heat to penetrate the load.

## CHEMICAL DISINFECTION

To disinfect linen chemically the process and products used must give results that are equal to that obtained by the thermal process.

The aerobic bacterial count must be less than 1 microorganism/cm<sup>2</sup>.

## THE SOLUTION

Cleveland VALIDATE Laundry Detergent.

VALIDATE complies with the Australian Standard for disinfecting in the wash.

VALIDATE will give the same result as thermal washing at a much lower temperature 40°C.

Ideal for washing personal clothing in smaller machines which cannot operate at higher temperatures.

## THE BENEFITS

Disinfects in the wash at lower temperatures

- Gives improved wash results.
- Removes offensive odours.
- Shorter wash times - gives labour cost saving
- Lower wash temperature - gives energy cost saving.

## NOW AVAILABLE IN PRE-MEASURED SOLUBLE SACHETS.

- No wastage
- No spills
- No handling chemical

---

CLEVELAND CLEANING SUPPLIES PTY. LTD. A.B.N. 61 001 175 748

Australia Wide: 1300 760 765

\*Head Office: 23 Forrester Street, Kingsgrove, N.S.W. 2208 Tel: (02) 9503 0100 Fax: (02) 9503 0199

Postal: P.O. Box 320, Kingsgrove, N.S.W. 1480 Australia

Sydney Tel: (02) 9503 0100 Melbourne Tel: (03) 9543 8998 Adelaide Tel: 1300 760 765

Canberra Tel: (02) 6280 7315 Dubbo Tel: (02) 6882 9433 Wagga Tel: (02) 6925 1303

Email: sales@cleveland.com.au Web: <http://www.cleveland.com.au>

